

CS-Pro MED

Therapeutic Volume Visualization

Five Versatile Standoffs

Enable precise treatment of a wide
range of indications

2mm

- Therapeutic Depth: 21mm
- Therapeutic Volume: 500mm³

5mm

- Therapeutic Depth: 25mm
- Therapeutic Volume: 750mm³

10mm

- Therapeutic Depth: 30mm
- Therapeutic Volume: 1000mm³

20mm

- Therapeutic Depth: 40mm
- Therapeutic Volume: 1500mm³

30mm

- Therapeutic Depth: 60mm
- Therapeutic Volume: 2000mm³

Note: The term "Therapeutic" defines the area
where the device produces at least 5 MPa.

